



TECHNICAL DATA SHEET

Medical face mask

(MNS) Mouth and nose protector made of lightweight non-woven fabric with extra soft ear loops and a shapeable nose clip for comfort.

KSI CARE SERIE

Medical face mask

KSi Article No.: 128-KSI-Care-Medizinische Gesichtsmaske

CERTIFICATIONS / PRODUCT-RELATED STANDARDS

BRI brand masks meet the following standards

Standards: Type IIR, CE-certified according to EN14683:2019

Type of individual packaging: 50 pieces / box

Type of EIS: Design board

Country of production: People's Republic of China

PRODUCT DESCRIPTION / MATERIAL

- Bacterial Filter Efficiency (BFE) \geq 98 %
- 3-layer filter
 - **Outer layer:** non-woven fabric
 - **Filter:** Melt blown material
 - **Inner layer:** non-woven fabric
- latex-free, fibreglass free, non-sterile
- Size: 9,5 x 17,5 cm

FEATURES

NR (non reusable) = For single use only

Breathable material allows easy inhalation and exhalation.

The mask has an optimal fit due to two elasticated elastic bands.

High wearing comfort through soft ear loops



AREAS OF APPLICATION

For medical and everyday use.

The mask covers mouth, nose and chin to act as a protective barrier to block the risk of direct transmission of infectious germs through coughing and sneezing.

PACKAGING

Quantity inner packaging	50 pieces per UVE
Quantity outer carton	2000 pieces in master carton
Full pallet	32,000
L x W x H Outer box in cm	55 x 41 x 43 cm
Volume in m ³ outer carton	0,097
Gross weight outer box in Kg.	9,5

PRODUCT FEATURES

Customs tariff number	63079098990
Minimum shelf life	2 years

SETUP INSTRUCTIONS



1 Before using the mask, wash hands thoroughly with soap or use disinfectant.



2 Place the mask on the nose and mouth with the coloured side facing outwards and the clip on top. Pull the rubber bands over the ears.



3 Press the nose clip to adjust it and have an optimal hold on the nose.



4 Adjusting the mask, make it cover the nose to the mandible, to get the best protective effect.



To take off the mask: Remove the mask from behind (do not touch the front); clean hands with alcohol hand rubs or soap and water.

Do you need more detailed information? We are happy to advise you!

The technical data are average production values and may vary slightly in individual cases. Subject to change without notice, no guarantee for printing errors or mistakes.